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SMALL Cartridge

4300-C2025 0.1μm pore size for the highest exchange rates. Activated fiber for attachment of matrix proteins, cytokines and antibodies. Ideal fiber for endothelial cell and hepatocyte culture.

MEDIUM Cartridges

4300-C2008 Low MWCO (5kd@50%) hydrophilic fiber for trapping smaller molecules. Suggested for recombinant proteins between 25kd and 100kd. Appropriate for suspension cell lines including CHO, HeLa and 293. Can support up to 10⁹ cells and produce 100μg/mL of recombinant protein in 15mL ECS.

4300-C2011 High MWCO (20kd@50%) hydrophilic fiber for trapping larger molecules, hybridoma culture and lymphocyte culture. 20kd MWCO allows TGF beta and TNF alpha to diffuse away while retaining antibodies. Can support up to 10⁹ cells and produce 5-50mg of monoclonal antibody every two days. Also used for medium scale adenovirus production at levels of 1- 5 x 10¹² pfu (plaque forming units).

LARGE Cartridges

4300-C2003 Low MWCO (5kd@50%) hydrophilic fiber for trapping smaller molecules. Suggested for recombinant proteins between 25kd - 100kd.

4300-C2018 High MWCO (20kd@50%) hydrophilic fiber for trapping larger molecules, hybridoma and lymphocyte cultures. 20kd MWCO allows TGF beta and TNF alpha to diffuse away while retaining antibodies. Can support up to 10¹¹ cells and produce 75-150mg of monoclonal antibody every two days. Appropriate for adherent suspension cell lines including CHO, HeLa and 293 cells. Can support up to 10¹¹ cells and produce 100μg/mL of recombinant protein in 70mL ECS.

X-LARGE Cartridges

35 Weaver Street, Scarsdale, New York 10583. U.S.A.

4300-C4005 The C4005 cartridge is intended for use in larger hollow fiber cell culture systems from other manufacturers and does not include a flow path stand or oxygenator tubing. Side ports have 3" of tubing capped with luer fittings, end ports are 3/8" hose barbs. High gross filtration rate and polysulfone fiber are superior to cellulose acetate for recombinant protein production.

4300-C4020 The C4020 cartridge is intended for use in larger hollow fiber cell culture systems from other manufacturers and does not include a flow path stand or oxygenator tubing. Side ports have 3" of tubing capped with luer fittings, end ports are 3/8" hose barbs with tubing and standard luer connectors on the ends. High gross filtration rate and polysulfone fiber are superior to cellulose acetate for recombinant protein and monoclonal antibody production.

What cell types have been cultured in the FiberCell Hollow Fiber bioreactor system.

Essentially if the cells can be grown in flask or other conventional systems then they will grow in a hollow fiber system. It is dependent upon your research and production goals as to whether the cells will behave in the desired fashion. The most common cells types successfully used are:

- Hybridoma cell lines of all species including NSO.
- HEK 293 both suspension and adherent
- CHO
- SP2
- HELA
- Hepatocytes
- Cancer cell lines
- COS cells
- HIV and B lymphocytes

35 Weaver Street, Scarsdale, New York 10583. U.S.A.

- And many others.

What are the advantages of hollow fiber cell culture?

Hollow fiber bioreactors support cells at 10-100X higher density than regular cell culture methods. This means the cells are in a more *in vivo* like environment and require less serum, can be more easily adapted to a serum free medium or can be supported with a simplified serum replacement like CDM HD.

Secreted products will be concentrated by the filter-like behavior of the hollow fibers, typically 100X higher concentration than that found with traditional bioreactors.

The effects of cytokines such as TGF-Beta or TNF Alpha can be controlled by the selection of the pore size of the fiber.

Hollow fiber bioreactors permit the culture and handling of large numbers of cells in a way that might not be practical using other methods in most laboratories.

Cells are bound to a porous support so they are free to grow in a post confluent fashion. Cells do not need to be split and can grow for extended periods of time. Hybridomas will typically produce antibody for 6 months or longer, CHO and 293 for 3-6 months of continuous production. The record is 2 years of continuous growth of a glioma cell line.

How does the pump work?

The FiberCell Duet pump uses a positive pressure displacement method that incorporates two one way check valves to drive the medium through the cartridge. This ensures reliable flow and long cartridge life.

How is gas controlled?

There is a loop of silicone tubing wrapped around the core of the cartridge flow path stand. Silicone tubing is very gas permeable and the gas composition of the medium will be the same as the gas composition in the incubator.

How is temperature controlled?

The Duet pump is designed to fit inside a standard CO2 incubator, the thin cord is designed to fit through the glass door.

What flow rate should I use?

The rate limiting factor in hollow fiber cell culture is the low partial pressure of oxygen due

35 Weaver Street, Scarsdale, New York 10583. U.S.A.

to lower solubility at 37 degrees. For this reason generally flow rate should be at the higher levels, between 26-30 on the Duet Control Box for the larger cartridges (C2003 and C2018) and between 22-26 for the medium sized cartridges. At the initiation of culture you will want to use a somewhat lower flow rate in order to allow cytokines to concentrate around the cells.

I want to produce a monoclonal antibody. What cartridge should I use and how much monoclonal antibody can I produce?

Cat #C2011 or C5011 should be used. This MWCO allows TGF Beta to diffuse away while retaining the produced antibody.

C2011 will support up to 1-2 X10⁹ cells, this is equal to a one-liter culture or more.

- Produce 5-50mg of antibody every 2 days, average is 20mgs per harvest
- Continue to produce antibody for up to 6 months of continuous culture
- Consumes about 1 liter of cell culture medium every two days. To reduce medium consumption harvest more cells out.
- A single mouse used for ascites fluid production will produce 10-20mgs total antibody; each harvest is equal to a single mouse.
- FiberCell cartridge has about 4-5 times the production capacity of the CellLine flask per harvest.
- Endotoxin burden is 1/10 th that of ascites fluid production.
- Cannot be re-used but can be stored and re-inoculated with the same cell line. product in the supernatant will be your antibody, simplifying your purification.
- Hybridomas grown in the FiberCell cartridge can be more easily adapted to serum free cell culture medium or adapted to as low as 2% FBS. When CDM HD is used the only

C5011 Will support up to 2-4 X10¹⁰ cells ; this is equal to a two liter culture.

- Produce 10-100 mgs of antibody every 2 days, average is 40mgs per harvest
- Consumes about 2 liters of cell culture medium every day.
- Scale up cartridge from C2011.
- Endotoxin burden is 1/10 th that of ascites fluid production.

I want to produce a recombinant protein. What cells can I use and what sort of production can I expect?

Hollow fiber bioreactors can be used with any cell type that will grow in flask or spinner culture. Stable transfectants should be used to take full advantage of the long term production potential offered. CHO, 293, HEP G2 and many other cell types have been used. Insect cell culture is not ideal due to the transient nature of the culture but constitutive

35 Weaver Street, Scarsdale, New York 10583. U.S.A.

expression in S2 cells have had excellent results at lower induction agent concentrations

Recombinant Protein Production

The molecular weight of the protein to be produced determines the fiber MWCO to use.

For proteins larger than 100kd use 4300-C2011 or 4300-C2018

For proteins of 20kd to 100kd use 4300-C2008 or 4300-C2003

Production will be typically 100X that of flasks with harvested product concentrations between 100 micrograms and 300 micrograms per ml per day.

4300-C2011 Has 2200cm² of surface area, equal to 12.5 T175 flasks. Because of the way cells grow when attached to a fiber the total number of cells will be equal to 50-60 T175 flasks.

- If the protein of interest can be trapped by the molecular weight cut-off of the fiber it will reach a concentration 100 times that of the same cell line grown in flask culture when collected from the ECS (extra-capillary space)
- Average productivity is around 100µg/ml (of ECS volume) or up to 1mg per day.
- Proteins greater than 100kd in molecular weight will be trapped by this fiber.
- If the protein of interest is too small to be trapped (like cytokines and cell growth factors) it will reach a concentration of 10 times that of the same cell line grown in flask culture when collected from the reservoir bottle.
- Harvest every day instead of every two days
- Medium consumption same as for hybridomas.
- Reduction in serum concentration or easier adaptation to serum-free medium makes purification of proteins or identification of low concentration growth factors easier.

4300-C2018 Has 1.2m² of surface area. This is equal to 68 T175 flasks. Because of the way that cells grow on hollow fibers this will be equal to over 400 T175 flasks.

- Medium consumption will be as high as 4 liters per day
- Can harvest 5-10 mg of protein per day

What medium should I use and should I use serum?

Any cell culture medium that is used in flasks can be used with a hollow fiber bioreactor. However, there are some special considerations to keep in mind.

- The high cell density allows the reduction of the amount of serum to 2% and can

35 Weaver Street, Scarsdale, New York 10583. U.S.A.

facilitate the adaptation to serum free mediums. This also means that cells can be supported using a simplified replacement for serum such as CDM HD offered by FiberCell Systems.

- Protein free mediums, such as CDM HD, can be used but keep in mind that protein free mediums generally do not contain any attachment factors. Generally we want the cells to attach to the fiber so it is preferred to seed the cells in the presence of serum and then adapt to serum free medium or CDM HD once the cells have reached high density (i.e. consuming 1 gram or more of glucose per day)
- Generally you want to avoid a medium such as RPMI due to the low (2.5 grams per liter) concentration of glucose. This simply means that you would need to change the medium more often, an inconvenience.
- Cells are growing in a stable, shear free environment. The use of surfactants such as pluronic F60 is not required. CDM HD does not contain any surfactants or other membrane protectants.

Should I use antibiotics?

Unless you have a compelling reason to not use antibiotics in your culture medium FiberCell Systems recommends the use of standard concentrations of antibiotics. We can fix anything, as long as you don't contaminate the cartridge. Antibiotics can help prevent the occasional lapse in technique from spoiling a culture. If your concern is endotoxin, and you wish for the absolute lowest levels of endotoxin then it is recommended that you work without antibiotics. Antibiotics can shield an infection but permit endotoxin to accumulate.

What other equipment do I need?

The only other piece of equipment that may not be part of a standard cell culture laboratory would be a way to measure glucose. There are some fancy and expensive machines to measure glucose out on the market but a simple glucometer like the ones used by diabetics and available at just about any drugstore or pharmacy will do the trick. Keep in mind that their readings won't be that accurate above 3.5 grams per liter so it is important to know the starting glucose concentration of your medium.

What are the advantages of hollow fiber cell culture for endothelial cells?

For endothelial cell culture under chronic shear use C2025. This cartridge allows extra-cellular matrix proteins to be attached to the fiber. One of the few ways to grow endothelial cells (cells that line the interior of blood vessels) under conditions of medium flowing over them (like they experience in the body). When grown under these conditions they behave in a much more in vivo like manner. They lay down flat and form a monolayer, form tight junctions, and certain genes are turned on in response to this shear stress that are not expressed in static culture. Other cell types such as vascular smooth

35 Weaver Street, Scarsdale, New York 10583. U.S.A.

muscle and brain glial cells can be co-cultivated with these endothelial cells.

C4300-2025 will yield about 100µg of total RNA for gene expression analysis.

What are the advantages of using the hollow fiber system for in vitro toxicology?

Hollow fiber bioreactor cartridges from FiberCell Systems offer a simple way to set up a two-compartment model for *in vitro* toxicology with higher levels of reproducible control to complex growth, infection, treatment, and sampling regimens. This system permits more realistic simulation of in vivo drug effects in a dynamically controlled system providing data that more accurately reflects biological responses. The design is fully disposable and will take into account the potential use of weaponized pathogens and genetically modified organisms

What organisms have been used with this system?

- Bacteria including tuberculosis, anthrax, plague, and MDR staph aureus to name a few.
- Viruses including HIV
- Tumor cell lines including breast cancer.

What drugs have been used?

Any type of therapeutic compound can be used. FiberCell offers two fiber types, polysulfone and cellulosic. The polysulfone is preferred because the flux rate and therefore equilibration of drug across the fiber is quite rapid and the geometry of the fibers results in even distribution of the fiber bundle inside the housing. Cellulosic fibers are generally used when the compound to be tested is highly non-polar which can result in significant non-specific binding of the compound to the fiber. Cellulosic fibers will have much lower non-specific binding but also lower flux rates so somewhat slower equilibration times across the fiber.

Advantages of the hollow fiber system:

- Closed bio-safe system
- Organism load can be high enough to match human infections. A high starting number is required to uncover the emergence of drug resistance.
- Drug pharmacokinetics can be exactly modeled on human profile
- Many experiments can be run simultaneously
- Complex systems such as two drug or two cell type cultures can be easily set up.

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1. Knazek RA, Gullino PM, Kohler PO, Dedrick RL. Cell culture on artificial capillaries: an approach to tissue growth in vitro. *Science* 1972; 178(56):65-6.
2. Gardner TA, KoSC, Yang L, Cadwell JJS, Chung LWK, Kao C. Serum-free recombinant production of adenovirus using a hollow fiber capillary system. *BioTechniques* 2001; 30(2):422-8.
3. Isayeva T, Kotova O, Krasnykh V, Kotov A. Advanced methods of adenovirus vector production for human gene therapy: roller bottles, microcarriers, and hollow fibers. *BioProcessing J* 2003; 2(4):75-81.
4. Li T, Glushakova S, Zimmerberg J. A new method for culturing *Plasmodium falciparum* shows replication at the highest erythrocyte densities. *J Infect Dis* 2003; 187:159-62.
5. Ott MJ, Ballermann BJ. Shear stress-conditioned, endothelial cell-seeded vascular grafts: improved cell adherence in response to in vitro shear stress. *Surgery* 1995; 117(3):334-9.
6. Waybill PN, Chinchilli VM, Ballermann BJ. Smooth muscle cell proliferation in response to co-culture with venous and arterial endothelial cells. *J Vasc Interv Radiol* 1997; 8(3):375-81.
7. Redmond E, Cahill P, Sitzman J. Flow mediated regulation of endothelin receptors in co-cultured vascular smooth muscle cells: an endothelium-dependent effect. *J Vasc Res* 1997:895-902.

New Developments in Hollow-Fiber Cell Culture

John J.S. Cadwell

Abstract: The growth of genetic engineering and hybridoma technologies has been a strong impetus for the investigation of new techniques for cell culture. The search for a simple, cost-effective cell culture system has resulted in the birth of a new area of biotechnology: the use of bioreactors to grow cells on a large scale. A bioreactor based on hollow-fiber bundles is shown to be an ideal method for culturing cells in any laboratory. Classical applications of hollow-fiber systems are described, in addition to some more recent ones.

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HOLLOW FIBER BIOREACTORS AND RELATED PRODUCTS

Designed for the following applications:


Monoclonal Antibody Production

Recombinant Protein Production

Endothelial Cell Culture


Lymphocyte Culture and Expansion

In Vitro Toxicology

FiberCell Systems Inc.

a better way to grow cells

www.fibercellsystems.com

What is the FiberCell® System advantage?

- Provides 10X the gross filtration rate of cellulosic hollow fibers.
- Improves cell viability.
- Eliminates dead spots with uniform fiber spacing and “wavy” fibers within the cartridge.
- Dialyzes waste and inhibitory cytokines away from cells.
- Stabilizes the PH of the cell culture medium.
- Enhances bio-safety with closed system.
- Handles large cell numbers (up to 10^{11}).
- Concentrates secreted proteins and antibodies up to 100X in the small volume of the extra-capillary space.
- Provides flexible control options with unique variable rate pump system.
- Ensures long life of the cartridge and facilitates nutrient and waste exchange across the fiber with positive-pressure displacement pumping system.


ABOVE: Cross-section of hollow fibers in FiberCell® Systems hollow fiber bio-reactor. High-density cell growth is shown in extra-capillary space.

- Cells grow on and around hollow fibers.
- Fiber geometry is optimized for both adherent and suspension cell types.
- Small molecules such as lactate and glucose can easily cross the fiber.
- Large molecules such as monoclonal antibodies and proteins are retained and concentrated in the small volume of the extra capillary space.

THE FIBERCELL® SYSTEM

The hollow fiber cartridge used in the FiberCell System has wavy fibers to eliminate dead space.

Enlarged cross section showing only three wavy fibers.

Cross section of cartridge showing the fiber ends potted in polyurethane so that medium flows only through the inside of the fiber.

THE FIBERCELL® SYSTEM: Hollow Fiber Cell Culture

A fundamentally different approach to cell culture is now available. Modeled after the mammalian circulatory system hollow fiber cell culture offers the most *in-vivo* like manner to grow cells in any laboratory.

Hollow fibers are small tube-like filters approximately 200 microns in diameter whose molecular weight cut-off can be between 5kd and .1 μ m. These fibers are sealed into a cartridge shell so that cell culture medium pumped through the end of the cartridge will flow through the inside of the fiber while the cells are grown on the outside of the fiber. These fibers then create a semi-permeable barrier of defined molecular weight cut-off (MWCO) between the compartment in which the cells are growing and the medium is flowing. Since the cells are attached to a porous support (the hollow fiber) rather than a non-porous plastic dish nutrients are delivered from the bottom layer of cells on upwards. Splitting of the cells is not required and cultures can be maintained for many months of continuous production. When the secreted protein is retained in the extra-capillary space it will accumulate to a concentration of up to 100 times higher than with conventional flask or roller bottle culture.

Hollow fibers also provide a tremendous amount of surface area in a small volume. Cells grow on and around the fibers at densities of greater than 1×10^8 per ml. Hollow fiber cell culture is the only means to culture cells at *in vivo* like cell densities. Cell culture at high densities has the following benefits:

- Reduces serum requirements and facilitates adaptation to serum-free medium, including CDM-HD
- Increases the concentration of secreted product by 10 to 100 times
- Viral and parasitic infections proceed rapidly

FiberCell® Systems Hollow Fiber cell culture cartridges are the method of choice for the production of 10mg to gram quantities of recombinant proteins and monoclonal antibodies and the culture of 10^9 to 10^{11} cells.

Protocols are available for monoclonal antibody production, recombinant protein production using CHO, 293 and other cell types, lymphocyte expansion, conditioned medium and endothelial cell culture under defined shear stress.

A single FiberCell® Systems hollow fiber cartridge is equivalent to harvesting 20-200 roller bottles daily.

SUGGESTED APPLICATIONS

MONOCLONAL ANTIBODY PRODUCTION

The FiberCell® Systems hollow fiber cell culture system is the method of choice for the production of 100mg to 2 grams of antibody per month. The 20kd MWCO fiber allows the inhibitory factor TGF beta to diffuse away from the cells while retaining the monoclonal antibody in the small volume of the extra-capillary space.

- C2011 will produce between 10mg and 50mg of antibody every two days
- Low harvest volume of 15 to 70 mls
- Cartridges can be maintained for over 6 months of continuous production
- Endotoxin burden is 10X lower than ascites fluid, ideal for animal injections
- Easy adaptation to reduced serum, CDM-HD, or serum free medium
- Chimeric, humanized and non-murine antibodies can be easily produced


SECRETED PROTEIN PRODUCTION

Current methods for the production of secreted biologicals such as recombinant proteins, cytokines and conditioned medium involve the use of inefficient *in-vitro* culture systems such as roller bottles or large numbers of tissue culture flasks. FiberCell® Systems hollow fiber culture systems facilitate the easy culture of large numbers of cells ($1-5 \times 10^{10}$) replacing roller bottles or hundreds of flasks while concentrating the secreted protein 100X or more.


- Produce 1-10 mgs of recombinant protein per day
- Secreted proteins can be 100 times more concentrated vs. tissue culture supernatant
- Reduced serum requirements facilitates purification
- Small harvest volume for easy handling
- No splitting of cells required, cartridge maintenance just 15 minutes a day
- Cultures can be maintained for several months of production
- Optimal cell culture conditions can result in improved protein assembly and folding.

LYMPHOCYTE CULTURE AND EXPANSION

FiberCell® Systems hollow fiber cell culture cartridges are ideal for the expansion of PBLs, lymphocyte subsets, HIV production and to perform large scale viral transductions of lymphocytes. The uniform distribution of the fiber bundle within the cartridge is optimum for suspension cell culture while the extremely high gross filtration rate of the hydrophilic polysulfone fiber insures adequate nutrient delivery and waste removal for rapidly growing cell types.

A) MW Standards
B, D, E, F, G) Harvested Antibody
C) Circulating serum free medium


Comparison of a purified recombinant human hexamerized IgG produced from CHO cells grown in either a 150cm² tissue culture flask (upper panel) or a FiberCell polysulphone cartridge (lower panel). Gel filtration chromatography reveals incomplete polymerization of the hexamerized IgG produced in the flask. Both a highly polymerized hexamerized IgG (peak A.) and non-polymerized (peak B.) are observed. In contrast, the same protein produced in a FiberCell cartridge is expressed almost entirely in highly polymerized form (peak C.) 478 mg of purified protein were produced in 2 months in a volume of less than 5 liters using catalog #2018.

Data courtesy of Dr. Jim Arthos - Bethesda, MD.

IN VITRO TOXICOLOGY

Hollow fiber cartridge systems from FiberCell Systems have been used to model pharmacologic bio-availability for different antimicrobial drugs and also to mimic dosage profiles that generate resistant organisms. Anti-cancer agents, anti-fungals, antibiotics, anti-virals (including anti-HIV) and anti-parasitic drugs have all been tested in the hollow fiber system. This system mimics *in vivo* drug dosage profiles in a rigorously controlled *in vitro* system composed of a hollow fiber cartridge and associated pumps and reservoirs. Bio-availability can be controlled in a precise fashion and any dosage profile can be evaluated.

ENDOTHELIAL CELL CULTURE

The FiberCell® Systems PS+ module (cat#4300-C2025) is especially designed for the culture and study of endothelial cells under flow. The PS+ fiber can easily be coated with extra-cellular matrix proteins (and/or cytokines and antibodies) to permit the attachment of endothelial cells to the interior wall of the fiber. The microprocessor controlled FiberCell® Systems Duet pump can be programmed to produce consistent and defined amounts of shear stress by regulating the flow of medium over the cells. This allows the study of endothelial cells in a more physiologic environment when compared to other methods. As a result, typical cell morphology is retained, a monolayer is formed and tight junctions are created. These culture conditions more closely mimic *in vivo* cell growth.


Figure 7 Bovine Aortic endothelial cells grown on a matrix of fibronectin.

The microporous nature of the fibers and the ability to control the extra-cellular matrix provides an ideal system for cellular co-cultivation with other cell types such as vascular smooth muscle or neuroglia.

- Inner surface area of the fibers is equal to a T75 flask
- Shear stress of .5 dynes/cm² to 25 dynes/cm² can be generated
- Approximately 100 micrograms of RNA can be isolated from each cartridge
- Most types of microscopy can be performed including immunohisto-chemical techniques

THE FIBERCELL® DUET™

FiberCell® Duet™ Specifications

Dimensions	9.5 x 16.5 x 8" (w x l x h)
Weight	9.5 lbs
Voltage	100, 120, 220 (50 or 60 hz)
Flow Rate	1-200mls, continuously variable

The FiberCell® Duet™ The FiberCell Duet pump is engineered for flexibility and reliability. It is designed to support either one or two FiberCell Hollow Fiber Bioreactor cartridges and flow paths with up to two separate 2 liter medium reservoirs. The FiberCell Duet pump can generate a flow rate of 1 ml/min up to 200 mls/min for precise control of flow rate for endothelial cell culture applications and high maximum flow rates for most efficient production when using the larger C2003 and C2018 cartridges. It utilizes a new low voltage brushless motor for long term, low heat operation and operation with 100 volt, 120 volt and 220 volt electrical input. The Duet pump utilizes our proprietary positive-pressure displacement pumping systems that generates high flow rates without wear on the pump tubing permitting continuous operation of the cartridges for several months to one year or more. The ergonomic, low-rise design provides areas to accommodate up to two 2 liter medium reservoir bottles.

- Occupies ½ shelf, 1/3 height of a standard CO² incubator
- Separate medium reservoirs eliminates cross contamination
- Low voltage power supply to minimize heat generation
- Versatile, can support small, medium and large cartridges (C2025, C2008, C2011, C2003, C2018, C5011)
- High flow rate for optimum culture using large cartridges
- Full Two year limited warranty

P3202 FiberCell Duet Pump, Universal Voltage, 100V, 120V, 220V

RESERVOIR CAP

CDM H-D

The FiberCell® reservoir cap is a re-usable, autoclavable cap designed to fit onto a standard cell culture media bottle. The 33mm cap (cat# A1005) will fit onto standard glass medium bottle and the 38mm cap (cat# A1006) will fit onto standard plastic medium bottles.

A1005	33mm Reservoir Cap Assy
A1006	38mm Reservoir Cap Assy

CDM-HD Serum Replacement is a chemically defined, protein free serum replacement that permits any basal medium such as DMEM to be used without serum. CDM-HD is designed specifically for the culture of cells at high density and optimized for use in our hollow fiber bioreactor systems. CDM-HD provides lot-to-lot consistency, simplifies purification and is an economical replacement for serum. It is available as a dry powder to make up one liter and is used at a concentration of 10%.

CDM-HD-1 Powder to make 1 Liter

FiberCell Systems offers additional technical support for our systems as well as for other hollow fiber bioreactor systems.

CARTRIDGE SPECIFICATIONS

Stock No.	Size	Surface Area	Fiber Type	Packing Density	ECS Volume	MWCO 50%	MWCO 95%	Maximum Cell#
C2025	Small	75cm ²	Activated PS	30%	1.5mL	0.1µm	0.1µm	10 ⁸
C2008	Medium	2100cm ²	low flux PS	50%	15mL	5kd	20kd	10 ⁹
C2011	Medium	2100cm ²	high flux PS	50%	15mL	20kd	100kd	10 ⁹
C2003	Large	1.2m ²	low flux PS	50%	70mL	5kd	20kd	5 × 10 ¹⁰
C2018	Large	1.2m ²	high flux PS	50%	70mL	20kd	100kd	5 × 10 ¹⁰
C3008	Medium	2000cm ²	cellulosic	38%	12mL	10kd	5kd	N/A
C4005	X-Large	2.5m ²	low flux PS	50%	150mL	5kd	20kd	10 ¹¹
C4020	X-Large	2.5m ²	high flux PS	50%	150mL	20kd	100kd	10 ¹¹
C5011	Medium	2100cm ²	high flux PS	50%	15mL	20kd	100kd	2 × 10 ⁹
Oxy 0.6	12mls	6000cm ²	Hydrophobic polyethylene	N/A	40mL	100kd	10kd	N/A

SUGGESTED APPLICATIONS FOR FIBERCELL® CARTRIDGES

C2025 0.1µm pore size for the highest exchange rates. Activated fiber for attachment of matrix proteins, cytokines and antibodies. Ideal fiber for endothelial cell and hepatocyte culture.

C2008 Low MWCO (5kd @ 50%) hydrophilic fiber for trapping smaller molecules. Suggested for recombinant proteins between 20kd and 100kd. Appropriate for suspension and adherent cell lines including CHO, HeLa and 293. Can support up to 10⁹ cells and produce 1-2 mgs of recombinant protein per day.

C2011 High MWCO (20kd @ 50%) hydrophilic fiber for trapping larger molecules, hybridoma culture and lymphocyte culture. 20kd MWCO allows TGF beta and TNF alpha to diffuse away while retaining antibody. Can support up to 10⁹ cells and produce 5-50 mg of monoclonal antibody every two days. Suggested for recombinant proteins larger than 100kd and will produce 1-2 mgs of protein per day.

C2003 Low MWCO (5kd @ 50%) hydrophilic fiber for trapping smaller molecules. Suggested for recombinant proteins between 20kd - 100kd. Appropriate for suspension cell lines including CHO, HeLa and 293 cells. Can support up to 5 × 10¹⁰ cells and produce 5-10 mgs of recombinant protein per day.

C2018 High MWCO (20kd @50%, 95% @ 100kd) for trapping larger molecules in the range of 100kd and larger such as recombinant proteins from CHO, HEK 293 and other cell types. Can support up to 5 × 10¹⁰ cells and produce 5-10 mgs of protein per day.

C3008 utilizes a cellulosic fiber with a nominal MWCO of 5KD. C3008 is intended primarily for use for in vitro toxicology applications where organic based drugs may exhibit non-specific binding to other fiber types.

C4005 The C4005 cartridge is intended for use in larger hollow fiber cell culture systems from other manufacturers and does not include a flow path stand or oxygenator tubing. Side ports have 3" of tubing capped with luer fittings, end ports are 3/8" hose barbs. High gross filtration rate and polysulfone fiber are superior to cellulose acetate for recombinant protein production.


C4020 The C4020 cartridge is intended for use in larger hollow fiber cell culture systems from other manufacturers and does not include a flow path stand or oxygenator tubing. Side ports have 3" of tubing capped with luer fittings, end ports are 3/8" hose barbs. High gross filtration rate and polysulfone fiber are superior to cellulose acetate for recombinant protein and monoclonal antibody production.

C5011 Same sized cartridge as the C2011 but optimized for maximum monoclonal antibody production by increasing flow rate and doubling the amount of oxygenation tubing


Oxy 0.6 The Oxy 0.6 is a replacement oxygenator cartridge for 1.0 m2 oxygenators available from other manufacturers. It utilizes a new fiber with better gas transfer and reduced medium weepage. The fibers are also uniformly spaced apart using an inert matting system further improving gas transfer. Module comes with luer connections on both the end ports and side ports. Due to the geometry of the fiber medium is flowed through the ECS and gas is sparged through the inside of the fibers.

C2025

C2008
C2011

C2003
C2018


OUR MISSION

At FiberCell® Systems Inc. we are committed to providing easy-to-use hollow fiber cell culture systems that allow researchers and biotechnology companies to produce quantities of cells and cell products not possible using traditional cell culture methods. We are constantly creating new protocols and methods for the use of our hollow fiber bioreactors so that any laboratory can take advantage of the benefits offered by FiberCell® Systems Inc.

John J. S. Cadwell
President and CEO

To order contact:

FiberCell Systems Inc.
a better way to grow cells

905 W. 7th Street #334
Frederick, Md. 21701
+86-400-6762467 028-86653817

www.ahmedical.com; E-mail: ahmedical@verizon.net